

BRILLIANT GREEN AGAR BASE, HUMAN

A selective and differentiation medium for the isolation of *Salmonella* spp. (including *S. typhi*) from clinical specimens.

Dehydrated media	
Code number:	500 g: BGH20500, 5 kg: BGH25000
Colour:	Beige
Appearance:	Homogeneous hygroscopic powder
pH before autoclaving (25 °C):	7,1 – 7,5

Direction: Suspend **21,5 g** in 500 ml of distilled water and heat with frequent agitation until the medium boils well. Sterilise by autoclaving at 121 °C for 15 minutes. Cool to 50 °C and add aseptically **10 drops (0,5 ml) Brilliant Green Solution, Sterile (BGS80030-1D)**. Mix well before pouring.

Prepared media	
Bottled media bases:	100 ml: BGH30100, 500 ml: BGH30500
Plated media:	55 mm: BGH50055, 90 mm: BGH50090
Colour of bottled media:	Pink
Colour of plated media:	Bluish
pH (25 °C):	7,2 – 7,4

Direction: Supplement the melted bottled media bases according to the direction of the dehydrated media and dispense aseptically into sterile Petri-dishes. Media in Petri-dishes are ready to use.

FORMULA in g/l

Peptones	16,50
Lactose	10,00
Sucrose	1,00
Glucose	0,50
Acid fuchsin	0,08
Agar	15,00

Note: The typical formula can be adjusted to obtain optimal performance.

Storage conditions: Store the dehydrated media tightly closed in a dry place at room temperature. Store the bottled media protected from light at room temperature. Store the plated media protected from light at 2-8 °C. Use before the expiry date on the label.

Quality control:

Test strains	Incubation temp: 37 °C	Growth	Incubation time: 48 h
<i>Salmonella typhi</i> ATCC 19430		Good, colourless colonies	
<i>Escherichia coli</i> ATCC 25922		Slightly inhibited, red colonies	
<i>Proteus mirabilis</i> ATCC 29906		Partially inhibited, colourless colonies without swarming	
<i>Enterococcus faecalis</i> ATCC 29212		Inhibited	

In vitro diagnostic – for professional use only!