

## BRILLIANT GREEN SOLUTION, STERILE

A sterile brilliant green solution for the isolation of *Salmonella* spp. (including *S. typhi*) from clinical specimens.

Description	
Code number:	30 ml: BGS80030-1D
Colour:	Dark green
Appearance:	Transparent solution

**Direction:** Suspend **21,5 g of Brilliant Green Agar Base, Human (BGH20500)** in 500 ml of distilled water and heat with frequent agitation until the medium boils well. Sterilise by autoclaving at 121 °C for 15 minutes. Cool to 50 °C and add aseptically **10 drops (0,5 ml) Brilliant Green Solution, Sterile (BGS80030-1D)**. Mix well before pouring.

### FORMULA

Brilliant green	0,1%
Solvent	Sterile distilled water

**Storage conditions:** Protected from light, at 2-8 °C.

**Warning!**

Use before the expiry date on the label.

**In vitro diagnostic – for professional use only!**