

COLUMBIA CNA AGAR BASE

A selective medium for the isolation and differentiation of Gram-positive micro-organisms.

Dehydrated media	
Code number:	500 g: CNA20500, 5 kg: CNA25000
Colour:	Yellowish
Appearance:	Homogeneous hygroscopic powder
pH before autoclaving (25 °C):	7,1 – 7,5

Direction: Suspend **42 g** in 950 ml of distilled water and heat with frequent agitation until the medium boils well. Sterilise by autoclaving at 115 °C for 15 minutes. Cool to 50 °C and add aseptically **50 ml of sterile defibrinated sheep blood**. Mix well before pouring.

Warning!

The medium is heat sensitive.
No further sterilisation is necessary or desirable.

Prepared media	
Bottled media bases:	100 ml: CNA30100, 500 ml: CNA30500
Plated media:	55 mm: CNA50055, 90 mm: CNA50090
Colour of bottled media bases:	Yellowish
Colour of plated media:	Ruby red
pH (25 °C):	7,2 – 7,4

Direction: Supplement the melted bottled media bases according to the direction of the dehydrated media and dispense aseptically into sterile Petri-dishes. Media in Petri-dishes are ready to use.

FORMULA in g/l

Nutrient substrate (peptones, extracts)	23,00
Starch, soluble	1,00
Sodium chloride	5,00
Nalidixic acid	0,01
Colistin	0,01
Agar	13,00

Note: The typical formula can be adjusted to obtain optimal performance.

Storage conditions: Store the dehydrated media tightly closed in a dry place at room temperature. Store the bottled media protected from light at room temperature. Store the plated media protected from light at 2-8 °C. Use before the expiry date on the label.

Quality control:

Test strains	Incubation temp: 37 °C	Growth	Incubation time: 24 h
<i>Streptococcus pneumoniae</i> ATCC 49619		Good, alpha haemolysis (under micro-aerobic conditions)	
<i>Streptococcus pyogenes</i> ATCC 19615		Good, beta haemolysis (under micro-aerobic conditions)	
<i>Enterococcus faecalis</i> ATCC 29212		Good, without haemolysis	
<i>Escherichia coli</i> ATCC 25922		Inhibited	

References: Ellner et al. (1966) Am. J. Clin. Pathol. 45: 502.

In vitro diagnostic – for professional use only!