

GC SELECTIVE SUPPLEMENT, VCNT

FOR 500 ml OF SELECTIVE THAYER-MARTIN AGAR, VCNT

Freeze-dried mixture for the isolation of pathogenic *Neisseria* spp. This supplement prevents the swarming of *Proteus* spp.

Description	
Code number:	4 ml: VCT80004
Colour:	White
Appearance:	Homogeneous lyophilisate

Direction: Suspend **19,5 g of GC Agar Base (GCA20500)** in 460 ml of distilled water and heat with frequent agitation until the medium boils well. Sterilise by autoclaving at 121 °C for 15 minutes. Cool to 50 °C and add aseptically **35 ml of sterile defibrinated blood** and “chocolate” by heating at 80 °C for 10 min. Cool to 50 °C. Dissolve the contents of **one vial of Growth Factor Mixture Hydration Fluid** with 5 ml of sterile distilled water and add aseptically to the **Growth Factor Mixture (GFM80005)**. Mix well and add aseptically to the medium. Dissolve the contents of **one vial of GC Selective Supplement, VCNT (VCT80004)** with 4 ml of sterile distilled water and add aseptically to the above. Mix well before pouring.

FORMULA [mg/vial]

Colistin	3,75
Trimethoprim	2,50
Vancomycin	1,50
Nystatin	1,50

Storage conditions: Protected from light, at 2-8 °C.

Warning!

Use before the expiry date on the label.

In vitro diagnostic – for professional use only!