

STAPH/STREP SELECTIVE SUPPLEMENT

FOR 500 ml OF COLUMBIA CNA AGAR
FOR 500 ml OF TODD-HEWITT SELECTIVE BROTH

Freeze-dried mixture for the isolation of *Staphylococcus* and *Streptococcus* spp.

Description	
Code number:	4 ml: SHS80004
Colour:	White
Appearance:	Homogeneous lyophilisate

Direction for Columbia CNA Agar: Suspend **21 g Columbia Blood Agar Base (COL20500)** in 475 ml of distilled water and heat with frequent agitation until the medium boils well. Sterilise by autoclaving at 121 °C for 15 minutes. Cool to 50 °C and add aseptically **25 ml of sterile defibrinated sheep blood** and the contents of **one vial of Staph/Strep Selective Supplement (SHS80004)** reconstituted with 4 ml of sterile distilled water. Mix well before pouring.

Direction for Todd-Hewitt Selective Broth: Suspend **18,5 g of Todd-Hewitt Broth (THB20500)** in 500 ml of distilled water and heat gently to dissolve the medium completely. Sterilise by autoclaving at 121 °C for 15 minutes. Cool to 50 °C and add aseptically the contents of **one vial of Staph/Strep Selective Supplement (SHS80004)** reconstituted with 4 ml of sterile distilled water. Mix well and dispense aseptically into sterile final containers.

FORMULA [mg/vial]

Nalidixic acid	5,0
Colistin	5,0

Storage conditions: Protected from light, at 2-8 °C.

Warning!

Use before the expiry date on the label.

In vitro diagnostic – for professional use only!