

BILE SALT No.3

Bile salt No.3 is prepared by refinement of bile salt in order to meet the demand for use as a selective agent.

Code number:	100 g: BBS10100, 500 g: BBS10500
Colour:	White
Appearance:	Fine powder, easily soluble in water

Physico-chemical characteristics

Parameter	Specification
Solubility in water 2%	Complete
pH (2% autoclaved solution)	7,5 – 9,0
Loss on drying	≤ 5 %
Residue on ignition	≤ 15 %
Sodium cholate	45 – 55 %
Sodium deoxycholate	45 – 55 %

Storage conditions: Protected from light, at room temperature.

Warning!

Hygroscopic product. Avoid heat and moisture.

In vitro diagnostic raw material – for professional use only!