

CM BROTH FOR AEROBIC ENRICHMENT

A highly nutritious medium for the cultivation of fastidious micro-organisms.

Prepared media	
Tubed media:	150 x 15 mm: CMA40010 (1,5 - 2 cm meat particles + 6 ml broth)
Colour:	Brown + yellowish
pH (25 °C):	7,0 - 7,2

Direction: Media in tubes are ready to use.

FORMULA in g/l

Lean meat	400,00
Peptones	35,00
L-Cysteine	0,40
Glucose	5,00
Sodium chloride	3,30
Vitamins	0,01
Buffers	5,00

Note: The typical formula can be adjusted to obtain optimal performance.

Storage conditions: Store the media protected from light at 2 - 8 °C.

Quality control:

Test strains	Incubation temp: 37 °C	Growth	Incubation time: 24 h
<i>Streptococcus pneumoniae</i> ATCC 49619		Good	
<i>Staphylococcus aureus</i> ATCC 29213		Good	

In vitro diagnostic – for professional use only!