

**CAMPYLOBACTER SELECTIVE SUPPLEMENT, BOLTON****FOR 500 ml OF BOLTON BROTH**

Freeze-dried mixture for the selective enrichment of *Campylobacter* spp.

<b>Description</b>	
Code number:	4 ml: CBS80004
Colour:	White
Appearance:	Homogeneous lyophilisate

**Direction:** Suspend **14 g of Bolton Broth Base (BOB20500)** in 470 ml of distilled water and heat gently to dissolve the medium completely. Sterilise by autoclaving at 121 °C for 15 minutes. Cool to 50 °C and add aseptically **25 ml of sterile lysed horse blood** and the contents of **one vial of Campylobacter Selective Supplement, Bolton (CBS80004)** reconstituted with 4 ml of sterile distilled water. Mix well and dispense aseptically into sterile test tubes.

**FORMULA [mg/vial]**

Cycloheximide	25
Cefoperazone	10
Trimethoprim	10
Vancomycin	10

**Storage conditions:** Protected from light, at 2-8 °C.

**Warning!**

Use before the expiry date on the label.

**In vitro diagnostic – for professional use only!**