

**DRBC PLUS SELECTIVE SUPPLEMENT****FOR 500 ml OF ROSE BENGAL DICHLORAN AGAR PLUS**

Freeze-dried mixture for the enumeration of yeasts and moulds.

<b>Description</b>	
Code number:	4 ml: RBS80004
Colour:	White
Appearance:	Homogeneous lyophilisate

**Direction:** Suspend **15,5 g Rose Bengal Dichloran Agar (RBD20500)** in 500 ml of distilled water and heat with frequent agitation until the medium boils well. Sterilise by autoclaving at 121 °C for 15 minutes. Cool to 50 °C and add aseptically the content of **one vial of DRBC Plus Selective Supplement (RBS80004)** reconstituted with 4 ml of sterile distilled water. Mix well before pouring.

**FORMULA [mg/vial]**

Chlortetracycline	25
Zinc sulphate heptahydrate	10
Copper sulphate pentahydrate	5

**Storage conditions:** Protected from light, at 2-8 °C.

**Warning!**

Use before the expiry date on the label.

**In vitro diagnostic - for professional use only!**