

## ETHANOL, 70%

Ethyl alcohol solution for examination of dairy products.

<b>Solutions:</b>	
Code number:	DID61000-70 - 1000 ml
Colour:	Water clear
Appearance:	Transparent, residue-free solution

**Procedure:** Use the reagent according to the regulation of relevant standard.

### FORMULA

Ethanol (containing 15% propyl alcohol)	70 %
Solvent	Distilled water

**Storage conditions:** Store the bottles protected from light at room temperature.

#### Warning!

Use the reagent only before the expiry date indicated on the label.

**In vitro diagnostics - for professional use only!**