

## FluoroBio® LSB

A selective enrichment medium for the detection of coliforms by a fluorogenic procedure.

<b>Dehydrated media</b>	
Code number:	500 g: LSM20500, 5 kg: LSM25000
Colour:	Yellowish
Appearance:	Homogeneous hygroscopic powder
pH before autoclaving (25 °C):	6,6 – 7,0

**Direction:** Suspend **35 g** in one litre of distilled water and heat gently to dissolve the medium completely. Dispense into test tubes fitted with Durham tube and sterilise by autoclaving at 115 °C for 15 minutes.

### Warning!

The medium is heat sensitive.  
No further sterilisation is necessary or desirable.

<b>Prepared media</b>	
Bottled media:	100 ml: LSM30100, 500 ml: LSM30500
Tubed media:	150 x 15 mm: LSM40010 (10 ml)
Colour:	Yellowish
pH (25 °C):	6,7 – 6,9

**Direction:** Dispense the bottled media aseptically into sterile test tubes fitted with Durham tube. Media in tubes are ready to use.

## FORMULA in g/l

Tryptose	19,4
Lactose	5,0
Sodium chloride	5,0
Sodium lauryl sulphate	0,1
MUG	0,1
Buffers	5,4

**Note:** The typical formula can be adjusted to obtain optimal performance.

**Storage conditions:** Store the dehydrated media tightly closed in a dry place at room temperature. Store the bottled and tubed media protected from light at room temperature. Use before the expiry date on the label.

### Quality control:

Test strains	Incubation temp: 37 °C	Growth	Incubation time: 24 h
<i>Escherichia coli</i>	ATCC 25922	Good, with gas production, fluorescence at 366 nm	
<i>Salmonella typhimurium</i>	ATCC 14028	Good, without gas production	
<i>Enterococcus faecalis</i>	ATCC 29212	Inhibited	

**References:** APHA (1976) Compendium of Methods for the Microbiological Examination of Foods  
Kilian and Bulow (1984) Acta Path. Micr. Scand. Sect. B. 84: 245.

**In vitro diagnostic – for professional use only!**