

# GSP SELECTIVE SUPPLEMENT

## FOR 500 ml OF GSP AGAR

Freeze-dried mixture for the detection and differentiation of *Pseudomonas* and *Aeromonas* spp.

<b>Description</b>	
Code number:	4 ml: GSU80004
Colour:	White
Appearance:	Homogeneous lyophilisate

**Direction:** Suspend **23 g of GSP Agar Base (GSP20500)** in 500 ml of distilled water and heat with frequent agitation until the medium boils well. Cool quickly to 50 °C and add aseptically the contents of one vial of **GSP Selective Supplement (GSU80004)** reconstituted with 4 ml of sterile distilled water. Mix well before pouring.

**Warning!**

The medium is heat sensitive.  
No further sterilisation is necessary or desirable.

### FORMULA [mg/vial]

Penicillin	70
------------	----

**Storage conditions:** Protected from light, at 2-8 °C.

**Warning!**

Use before the expiry date on the label.

**In vitro diagnostic – for professional use only!**