

ITC SELECTIVE SUPPLEMENT

FOR 500 ml OF YERSINIA (ITC) BROTH

Freeze-dried mixture for the detection of *Yersinia enterocolitica*.

Description	
Code number:	4 ml: YIS80004
Colour:	White
Appearance:	Homogeneous lyophilisate

Direction: Fill up **29 ml of DIASALM-MSRV Magnesium Chloride Solution (DSM80500)** to 500 ml with distilled water. Suspend **8 g of Yersinia (ITC) Broth Base (YIT20500)** and heat gently to dissolve the medium completely. Sterilise by autoclaving at 121 °C for 15 minutes. Cool to 50 °C and add aseptically the contents of **one vial of ITC Potassium Chlorate Solution, Sterile (YIP80005)** and **one vial of ITC Selective Supplement (YIS80004)** reconstituted with 4 ml of sterile distilled water. Mix well and dispensed aseptically into sterile test tubes.

FORMULA [mg/vial]

Irgasan	0,5
Ticarcillin	0,5

Storage conditions: Protected from light, at 2-8 °C.

Warning!

Use before the expiry date on the label.

In vitro diagnostic – for professional use only!