

LISTERIA SELECTIVE SUPPLEMENT, UVM II

FOR 500 ml OF LISTERIA SELECTIVE BROTH, UVM II

Freeze-dried mixture for the isolation of *Listeria monocytogenes*.

Description	
Code number:	4 ml: LU280004
Colour:	White
Appearance:	Homogeneous lyophilisate

Direction: Suspend **27,5 g of Listeria Enrichment Broth Base, UVM-Fraser (LEF20500)** in 500 ml of distilled water and heat gently to dissolve the medium completely. Sterilise by autoclaving at 121 °C for 15 minutes. Cool to 50 °C and add aseptically the contents of **one vial of Listeria Selective Supplement, UVM II (LU280004)** reconstituted with 4 ml of 1:1 mixture of ethanol and sterile distilled water. Mix well. Dispense aseptically into sterile final containers.

FORMULA [mg/vial]

Nalidixic acid	12,5
Acriflavine	10,0

Storage conditions: Protected from light, at 2-8 °C.

Warning!

Use before the expiry date on the label.

In vitro diagnostic – for professional use only!