

## M17 BROTH

A selective medium for the cultivation and enumeration of lactic streptococci.

Dehydrated media	
Code number:	500 g: M1B20500, 5 kg: M1B25000
Colour:	Yellowish
Appearance:	Homogeneous hygroscopic powder
pH before autoclaving (25 °C):	6,7 – 7,1

**Direction:** Suspend **42 g** in one litre of distilled water and heat gently to dissolve the medium completely. Dispense into final containers and sterilise by autoclaving at 121 °C for 15 minutes.

Prepared media	
Bottled media:	100 ml: M1B30100, 500 ml: M1B30500
Tubed media:	150 x 15 mm: M1B40010 (10 ml)
Colour:	Yellowish
pH (25 °C):	6,8 – 7,0

**Direction:** Dispense the bottled media aseptically into sterile final containers. Media in tubes are ready to use.

### FORMULA in g/l

Peptones	17,25
Lactose	5,00
Magnesium sulphate	0,25
Ascorbic acid	0,50
Sodium glycerophosphate	19,00

**Note:** The typical formula can be adjusted to obtain optimal performance.

**Storage conditions:** Store the dehydrated media tightly closed in a dry place at room temperature. Store the bottled and tubed media protected from light at room temperature. Use before the expiry date on the label.

#### Quality control:

Test strains	Incubation temp: 30 °C	Growth	Incubation time: 48 h
<i>Streptococcus lactis</i> ATCC 19435		Good	

**References:** Terzaghi and Sandine (1975) Applied Microbiol. 29: 807.

**In vitro diagnostic – for professional use only!**