

MALACHITE GREEN SOLUTION, STERILE

FOR 10 OR 30 LITRE OF MALACHITE GREEN BROTH

Malachite green solution for the preparation of Malachite Green Broth.

Description	
Code number:	10 ml: MSO80010, 30 ml: MSO80030
Colour:	Green
Appearance:	Transparent solution

Direction: Suspend **4,2 g Malachite Green Broth Base (MIB20500)** in 500 ml of distilled water and heat gently to dissolve the medium completely. Sterilise by autoclaving at 121 °C for 15 minutes. Cool to 50 °C and add aseptically **10 drops (0,5 ml) of Malachite Green Solution, Sterile (MSO80030)**. Mix well and dispense aseptically into sterile final containers.

FORMULA

Malachite green	1 %
Solvent	Sterile distilled water

Storage conditions: Protected from light, at room temperature.

Warning!

Use before the expiry date on the label.

In vitro diagnostic – for professional use only!