

**OGYE SELECTIVE SUPPLEMENT****FOR 500 ml OF OXYTETRACYCLINE GLUCOSE YEAST EXTRACT  
AGAR**

Freeze-dried mixture for the isolation of yeasts and moulds.

Description	
Code number:	4 ml: OGS80004
Colour:	White
Appearance:	Homogeneous lyophilisate

**Direction:** Suspend **19 g Oxytetracycline Glucose Yeast Extract Agar Base (OGY20500)** in 500 ml of distilled water and heat with frequent agitation until the medium boils well. Sterilise by autoclaving at 121 °C for 15 minutes. Cool to 50 °C and add aseptically the contents of **one vial of OGYE Selective Supplement (OGS80004)** reconstituted with 4 ml of sterile distilled water. Mix well before pouring.

**FORMULA [mg/vial]**

Oxytetracycline	50
-----------------	----

**Storage conditions:** Protected from light, at 2-8 °C.

**Warning!**

Use before the expiry date on the label.

**In vitro diagnostic - for professional use only!**