

PEPTONE WATER, DOUBLE BUFFERED

A pre-enrichment medium for isolation of *Salmonella* spp.

Dehydrated media	
Code number:	500 g: PWD20500, 5 kg: PWD25000
Colour:	Yellowish
Appearance:	Homogeneous hygroscopic powder
pH before autoclaving (25 °C):	7,0 – 7,4

Direction: Suspend **25 g** in one litre of distilled water and heat gently to dissolve the medium completely. Dispense into final containers and sterilise by autoclaving at 121 °C for 15 minutes.

Prepared media	
Bottled media:	100 ml: PWD30100, 500 ml: PWD30500
Tubed media:	150 x 15 mm: PWD40010 (10 ml)
Colour:	Yellowish
pH (25 °C):	7,1 – 7,3

Direction: Dispense the bottled media aseptically into sterile final containers. Media in tubes are ready to use.

FORMULA in g/l

Peptones	10
Sodium chloride	5
Potassium phosphate, monobasic	3
Sodium phosphate, dibasic (Equivalent to 18.0 g of Sodium phosphate dibasic dodecahydrate)	7

Note: The typical formula can be adjusted to obtain optimal performance.

Storage conditions: Store the dehydrated media tightly closed in a dry place at room temperature. Store the bottled and tubed media protected from light at room temperature. Use before the expiry date on the label.

Quality control:

Test strains	Incubation temp: 37 °C	Growth	Incubation time: 24 h
<i>Salmonella typhimurium</i> ATCC 14028		Good	

References: ISO 6579

In vitro diagnostic – for professional use only!