

PERFRINGENS SELECTIVE SUPPLEMENT, OPSP, A + B

FOR 500 ml OF PERFRINGENS OPSP AGAR

Freeze-dried mixture for the isolation of *Clostridium perfringens*.

Description	
Code number:	2x4 ml: POS80004
Colour:	White
Appearance:	Homogeneous lyophilisate

Direction: Suspend **23,5 g of Perfringens (OPSP) Agar Base (POB20500)** in 500 ml of distilled water and heat with frequent agitation until the medium boils well. Sterilise by autoclaving at 121 °C for 15 minutes. Cool to 50 °C and add aseptically the contents of **one vial each of Perfringens Selective Supplements, OPSP, A + B (POS80004)** both reconstituted with 4 ml of sterile distilled water. Mix well before pouring.

FORMULA [mg/vial]

Perfringens Selective Supplement, OPSP, A	
Polymyxin B	0,80
Oleandomycin	0,25
Perfringens Selective Supplement, OPSP, B	
Sulfadiazine	50,00

Storage conditions: Protected from light, at 2-8 °C.

Warning!

Use before the expiry date on the label.

In vitro diagnostic – for professional use only!