

PHYSIOLOGICAL SALT SOLUTION

A sterile physiological salt solution for microbiological purposes.

Prepared solutions	
Bottled solutions:	50 ml: PSO30050 100 ml: PSO30100 500 ml: PSO30500
Tubed solutions:	150 x 15 mm: PSO40010 (10 ml)
	150 x 15 mm: PSO40009 (9 ml)
	100 x 15 mm: PSO40005 (5 ml)
	100 x 12 mm: PSO40003 (3 ml)
	100 x 12 mm: PSO40002 (2 ml)
	100 x 12 mm: PSO40001 (1 ml)
	100 x 12 mm: PSO400005 (0,5 ml)
Ampouled solutions:	10ml ampoule: PSO40004-AI (4 ml)
Colour:	Water clear
pH (25 °C):	6,9 – 7,1

Direction: Dispense the bottled solutions aseptically into sterile final containers.

FORMULA in g/l

Sodium chloride	9
-----------------	---

Storage conditions: Store the bottled, tubed and ampouled solutions protected from light at room temperature. Use before the expiry date on the label.

In vitro diagnostic - for professional use only!