

**SHIGELLA SELECTIVE SUPPLEMENT****FOR 500 ml OF SHIGELLA SELECTIVE BROTH**

Freeze-dried mixture for the selective enrichment of *Shigella* spp.

<b>Description</b>	
Code number:	4 ml: SBS80004
Colour:	White
Appearance:	Homogeneous lyophilisate

**Direction:** Suspend **15 g Shigella Broth Base (SHB20500)** in 500 ml of distilled water and heat gently to dissolve the medium completely. Sterilise by autoclaving at 121 °C for 15 minutes. Cool to 50 °C and add aseptically the contents of **one vial of Shigella Selective Supplement (SBS80004)** reconstituted with 4 ml of sterile distilled water. Mix well and dispense aseptically into sterile test tubes.

**FORMULA [mg/vial]**

Cefsulodin	20,0
Novobiocin	0,3

**Storage conditions:** Protected from light, at 2-8 °C.

**Warning!**

Use before the expiry date on the label.

**In vitro diagnostic – for professional use only!**