

## STERILE DISTILLED WATER

A diluent for miniaturized microplate method for the multi-tube (MPN) procedure according to the EN/ISO 9308-3 (MUG/EC) and EN ISO 7899-1 (MUD/SF).

Prepared media	
Bottled media:	100 ml: STW30100, 500 ml: STW30500
Tubed media:	150 x 15 mm: STW40018 (18 ml)
Colour:	Water-clear
pH (25 °C):	6,8- 7,2

**Direction:** The diluents are ready to use.

**Storage conditions:** Store protected from light at room temperature. Use before the expiry date on the label.

**Quality control:** according to ISO 9308-3 (MUG/EC) and ISO 7899-1 (MUD/SF).

**References:** ISO 9308-3:2000; ISO 7899-1:2000

**In vitro diagnostic – for professional use only!**