

# TRITON X-100 SUPPLEMENT

## FOR 100 OR 500 LITRE OF A-1 BROTH

TRITON X-100 (polyethylene glycol tert-octylphenyl ether) for preparation of A-1 Broth.

Description	
Code number:	100 ml: TXS80100, 500 ml: TXS80500
Colour:	Colourless or slightly yellowish
Appearance:	Transparent liquid

**Direction:** Suspend **31 g of A-1 Broth Base (A1B20500)** in one litre of distilled water. Add **1 ml of Triton X-100 Supplement (TXS80100)**. Mix well and heat gently to dissolve the medium completely. Dispense into test tubes fitted with Durham tube and sterilise by autoclaving at 121 °C for 15 minutes.

### FORMULA

TRITON X-100 (bacteriological grade)	TXS80100 - 100 ml
TRITON X-100 (bacteriological grade)	TXS80500 - 500 ml

**Storage conditions:** Protected from light, at room temperature.

#### Warning!

Use before the expiry date on the label.

**In vitro diagnostic – for professional use only!**