

TWEEN 80 SUPPLEMENT

TWEEN 80 (Polyoxyethylene sorbitan monooleate; Polysorbate 80) for the preparation of some media.

Description	
Code number:	100 ml: TWS80100, 500 ml: TWS80500
Colour:	Colourless or slightly yellowish
Appearance:	Thick, transparent liquid. A few precipitation may occur but not detrimental.

Direction: Various. See the product information of the relevant media bases.

FORMULA

TWS80100	TWEEN 80 (bacteriological grade) - 100 ml
TWS80500	TWEEN 80 (bacteriological grade) - 500 ml

Storage conditions: Protected from light, at room temperature.

Warning!

Use before the expiry date on the label.

In vitro diagnostic – for professional use only!