

TETRATHIONATE IODINE-IODIDE SELECTIVE SUPPLEMENT

FOR 500 ml OF TETRATHIONATE BROTHS

Dehydrated mixture for the preparation of tetrathionate broths.

Description	
Code number:	10 ml: TTS80010
Colour:	Dark gray - black
Appearance:	Shiny pellets

Direction: Suspend the appropriate quantity of one of the tetrathionate broth bases in 500 ml of distilled water and heat gently to dissolve the medium completely. Cool to 50 °C and add aseptically the contents of **one vial of Tetrathionate Iodine-Iodide Selective Supplement (TTS80010)** reconstituted with 10 ml of sterile distilled water. Mix well and dispense aseptically into sterile test tubes.

Warning!

The medium is heat sensitive.

Do not heat after the addition of the supplement.

It is recommended to use the complete medium on the day of preparation.

FORMULA [g/vial]

Iodine	3,0
Potassium iodide	2,5

Storage conditions: Protected from light, at 2-8 °C.

Warning!

Use before the expiry date on the label.

In vitro diagnostic – for professional use only!