

CRYOTUBE

The vials, which are manufactured from high quality polypropylene, are suitable for storage of biological specimens at extremely low temperature.

Name:	Cryotube with hub and gasket, 5 ml, sterile
	for laboratory use only
Code number:	TUB93002

Specification:	
Packaging:	200 pcs/bag 3200 pcs/case
Packaging material:	PE bag
Colour:	Slightly opalescent
Volume:	5 ml
Material:	Special grade polypropylene
Heat resistance:	Autoclavable
Sterility:	Gamma irradiated
SAL:	10 ⁻³
Note:	-

Storage conditions: Store in a dry place at room temperature.

Quality control: The conformity of the parameters above is checked by the Quality Control Laboratory.