

CULTURE MEDIA TUBE, GLASS, WITH METAL SCREW CAP

Glass, disposable, flat bottomed tube with metal screw cap for culture media.

Name:	Culture Media Tube, Glass, with Metal Screw Cap
	for laboratory use only
Code number:	GTU91013

Specification:	
Colour of tube:	Water clear
Colour of cap:	Gold
Material of tube:	Glass
Material of cap:	Aluminium
Seal:	Silicon rubber
Sizes:	115 x 16 mm
Sterility:	Non-sterile

Storage conditions: Store in a dry place at room temperature.

Quality control: The conformity of the parameters above is checked by the Quality Control Laboratory.