

## DURHAM TUBE

Glass tube for the investigation of gas production of micro-organisms.

<b>Name:</b>	<b>Durham tube, glass</b>
	for laboratory use only
Code number:	TUB90004

<b>Specification:</b>	
Packaging:	100 pcs/bag
Packaging material:	PE bag
Colour:	Water clear
Material:	Glass
Heat resistance:	Autoclavable
Sterility:	Non-sterile
Note:	-

**Storage conditions:** Store in a dry place at room temperature.

**Quality control:** The conformity of the parameters above is checked by the Quality Control Laboratory.