

EMPTY RACK FOR TIPS

Empty rack for storage and sterilisation of pipette tips.

Name:	Empty rack for 96 pcs of 200 µl pipette tips.
	for laboratory use only
Code number:	PIT92801

Specification:	
Packaging:	100 pcs/case
Packaging material:	PE bag
Colour:	Yellow
Material:	Polypropylene
Heat resistance:	Autoclavable
Sterility:	Non-sterile
Note:	-

Storage conditions: Store in a dry place at room temperature.

Quality control: The conformity of the parameters above is checked by the Quality Control Laboratory.