

## FAECAL CONTAINER

Plastic, disposable, conically bottomed container, with skirt and with screw cap for sampling, transport and storage of faeces samples.

<b>Name:</b>	<b>Faecal container with red screw cap, conical bottom with skirt, 30 ml, labelled, sterile</b>
	for laboratory use only
<b>Code number:</b>	FEC90019

<b>Specification:</b>	
Packaging:	500 pcs/case
Colour of container:	Slightly opalescent
Colour of cap:	Red
Volume:	30 ml
Material:	High quality polypropylene
Sterility:	Gamma irradiated
SAL:	10 <sup>-3</sup>
Note:	-

**Storage conditions:** Store in a dry place at room temperature.

**Quality control:** The conformity of the parameters above is checked by the Quality Control Laboratory.