

## INOCULATING LOOP

Sterile loop with needle for the inoculation of micro-organism.

<b>Name:</b>	<b>Inoculating loop, 10 µl, sterile</b>
	for laboratory use only
<b>Code number:</b>	SLO90002

<b>Specification:</b>	
<b>Packaging:</b>	20 pcs/blister 1.000 pcs/unite 10.000 pcs/case
<b>Packaging material:</b>	Paper/foil blister
<b>Colour:</b>	Yellow
<b>Length:</b>	200 mm
<b>Loop volume:</b>	10 µl
<b>Material:</b>	Polystyrene
<b>Sterility:</b>	Irradiated
<b>SAL:</b>	10 <sup>-3</sup>
<b>Note:</b>	-

**Storage conditions:** Store in a dry place at room temperature.

**Quality control:** The conformity of the parameters above is checked by the Quality Control Laboratory.