

PASTEUR PIPETTE, 1 ml, STERILE

The plastic Pasteur pipette is the more convenient alternative of the glass Pasteur pipettes. The product is disposable, chemically inert, non-toxic and significantly reduces the possibility of the contamination between the samples.

Name:	Pasteur pipette, 1 ml, sterile, individually wrapped
	for laboratory use only
Code number:	PPI92001

Specification:	
Packaging:	100 pcs/unit 2000 pcs/case
Material:	LDPE (low density polyethylene)
Colour:	Translucent
Nominal fill volume:	≈ 1 ml
Volume of drops:	≈ 50 µl
Sterility:	Gamma irradiated
SAL:	10 ⁻³
Note:	-

Storage conditions: Store in a dry place protected from light at room temperature.

Quality control: The conformity of the parameters above is checked by the Quality Control Laboratory.