

## PASTEUR PIPETTE, 3 ml, NON-STERILE

The plastic Pasteur pipette is the more convenient alternative of the glass Pasteur pipettes. The product is disposable, chemically inert, non-toxic and significantly reduces the possibility of the contamination between the samples.

<b>Name:</b>	<b>Pasteur pipette, 3 ml</b>
	for laboratory use only
Code number:	PPI91002

<b>Specification:</b>	
Packaging:	500 pcs/unit 5000 pcs/case
Material:	LDPE (low density polyethylene)
Colour:	Translucent
Nominal fill volume:	≈ 3 ml
Volume of drops:	≈ 50 µl
Sterility:	Non-sterile
Note:	-

**Storage conditions:** Store in a dry place protected from light at room temperature.

**Quality control:** The conformity of the parameters above is checked by the Quality Control Laboratory.