

PIPETTE TIP

Micropipette tips, made from high quality polypropylene. Tips are manufactured with the highest precision to ensure accurate and repeatable results in dispensing.

Name:	Pipette tip 1000 µl (Gilson type)
	for laboratory use only
Code number:	PIT91805

Specification:	
Packaging:	500 pcs/bag 10.000 pcs/case
Packaging material:	PE bag
Colour:	Blue
Material:	Polypropylene
Heat resistance:	Autoclavable
Suitable pipettes:	Gilson, Nichiryo, Finnpiquette, Biohit etc.
Sterility:	Non-sterile
Note:	-

Storage conditions: Store in a dry place at room temperature.

Quality control: The conformity of the parameters above is checked by the Quality Control Laboratory.