

PIPETTE TIP

Micropipette tips, made from high quality polypropylene. Tips are manufactured with the highest precision to ensure accurate and repeatable results in dispensing.

Name:	Pipette tip 10 µl (Eppendorf type)
	for laboratory use only
Code number:	PIT91807

Specification:	
Packaging:	1.000 pcs/bag 50.000 pcs/case
Packaging material:	PE bag
Colour:	Colourless
Material:	Polypropylene
Heat resistance:	Autoclavable
Suitable pipettes:	Eppendorf, Socorex, Brand etc.
Sterility:	Non-sterile
Note:	-

Storage conditions: Store in a dry place at room temperature.

Quality control: The conformity of the parameters above is checked by the Quality Control Laboratory.