

SCREW CAP TUBE

Plastic, disposable, conically bottomed tube, with skirt and with screw cap for centrifuge, transport and storage of materials.

Name:	Screw cap tube, conical bottom, 50 ml, with skirt, sterile
	for laboratory use only
Code number:	TUB92161

Specification:	
Packaging:	25 pcs/bag 500 pcs/case
Packaging material:	PE bag
Colour of tube:	Slightly opalescent
Colour of cap:	Orange
Colour of printed graduation:	White
Volume:	50 ml
Material:	High quality polypropylene
Sterility:	Gamma irradiated
SAL:	10 ⁻³
Note:	-

Storage conditions: Store in a dry place at room temperature.

Quality control: The conformity of the parameters above is checked by the Quality Control Laboratory.