

SPECIAL CAP FOR CULTURE MEDIA TUBE

Flexible cap for closure of Ø15 mm culture media tubes.

Name:	Special cap for Ø 15 mm tubes, bulk
	for laboratory use only
Code number:	GTU91008

Specification:	
Packaging:	According to the customer request
Packaging material:	PE bag
Colour:	Blue
Material:	Special elastomer
Sterility:	Non-sterile
Note:	-

Storage conditions: Store in a dry place at room temperature.

Quality control: The conformity of the parameters above is checked by the Quality Control Laboratory.