

SWAB APPLICATOR

Sterile dry swab in tube with cap, for taking different samples.

| | |
|---------------------|---|
| Name: | Swab applicator in tube, wooden, sterile |
| | for laboratory use only |
| Code number: | SWA90001 |

| | |
|---------------------------|-----------------------------|
| Specification: | |
| Packaging: | 20 pcs/bag 2500 pcs/case |
| Packaging material: | PE bag with slide fastener |
| Material of tube and cap: | Polypropylene homopolymer |
| Material of stick: | Wooden |
| Dimensions of stick: | 150 x 2,5 mm |
| Material of tip: | Cotton |
| Diameter of tip: | ≈ 5 mm |
| Sterility: | Gamma irradiated |
| SAL: | 10 ⁻³ |
| Note: | - |

Storage conditions: Store in a dry place at room temperature.

Quality control: The conformity of the parameters above is checked by the Quality Control Laboratory.